

Shoulder Immobilizers

Secure Immobilization | Adjustable Comfort | Universal Fit



Product Overview

The Shoulder Immobilizer with Foam Straps provides secure and comfortable support during shoulder injury recovery. Its durable cotton sling, adjustable straps, and waist strap help maintain proper arm positioning while limiting shoulder rotation. The universal design fits either arm, and the included shoulder pad enhances comfort during extended use.

Features and Benefits

- Secure Support: Helps immobilize the shoulder throughout recovery.
- Adjustable Fit: Hook-and-loop closures on all straps allow personalized positioning.
- Universal Design: Suitable for either the right or left arm.
- Rotation Control: Waist strap helps prevent unwanted shoulder movement.
- Enhanced Comfort: Included shoulder pad improves patient comfort during extended use.

Product Specifications

- Product Type: Shoulder immobilizer
- Clinical Applications: Post-trauma shoulder subluxations and repetitive strain injuries
- Sling Material: Cotton
- Loop Material: Polypropylene
- Strap Material: Polyester
- Closure Type: Hook-and-loop
- Orientation: Universal—suitable for either arm

Packaging Details

Item No.	QTY	Size
1000-00	40/CS	Extra Small
1000-01	40/CS	Small
1000-02	40/CS	Medium
1000-03	40/CS	Large
1000-04	40/CS	Extra Large

