

Compression Knee Immobilizers



Secure Immobilization | Even Compression | Universal Fit

Product Overview

The Compression Knee Immobilizer provides secure stabilization and dependable support throughout knee injury recovery. Its contoured, reinforced stays deliver effective rigidity, while the soft foam-laminate construction enhances patient comfort. Adjustable elastic straps provide even compression and help maintain a secure, customized fit during the healing process.

Features and Benefits

- Rigid Support: Contoured, reinforced stays with a center “I” beam provide rigidity and a custom fit.
- Comfortable Construction: Soft foam-laminate material adds comfort and resists staining.
- Adjustable Sizing: The material can be easily trimmed to accommodate individual requirements.
- Even Compression: Adjustable elastic straps provide secure and balanced compression.
- Universal Design: Designed to fit most patients.

Product Specifications

- Product Type: Compression knee immobilizer
- Intended Use: Post-surgical immobilization, patellar fractures, and patellar dislocations
- Soft-Goods Material: Open-cell polyurethane foam with polyester-spandex jersey knit fabric
- Stay Material: Aluminum
- Strap Material: Polyester and spandex
- Support Design: Contoured, reinforced stays with a center “I” beam
- Fit: Universal

Packaging Details

Item No.	QTY	Size
2000-00	12/CS	12 in
2000-02	12/CS	16 in
2000-03	12/CS	20 in
2000-04	12/CS	24 in

