

Wrist Brace

Custom Comfort

Our patented Wrist Brace is the first off-the-shelf, custom fitting wrist support. Offers optimal support and immobilization while allowing unrestricted motion of the hand and thumb. Automatically adjusts and conforms to each patient's anatomy, allowing one size to fit nearly all patients.

Key Features:

- Standard size fits virtually all patients, reducing need for excess inventory.
- Panel automatically expands, providing a custom off-the-shelf fit.
- Auto-lock strap secures custom sizing.
- Cut low at the palmar crease to allow for 90-degree MP flexion.

MATERIALS USED IN PRODUCT:

Component	Material	
Main Body Softgoods	Polyester	
Stay Pocket	Nylon	
Palm Alum. Stay	Aluminum	
Side velcro fur	Nylon	
Doral Stay Pocket	Nylon Hook and Loop	
Dorsal Alu. Stay	Aluminum	
Ulnar/radial stays	Polypropylene	
Lace connection	Nylon	
Thumb velcro hook	Nylon	
Rivet	Steel	
Tape	Nylon	
Stockinet	Polyester	
Thread	Polyester	

WRIST BRACE

Product	Item no.	QTY.	Size
Wrist Brace	3000-11	50/CS	XS-LT
	3000-12	50/CS	XS-RT
	3000-10	50/CS	Standard-LT
	3000-20	50/CS	Standard-RT
	3000-21	50/CS	XL-LT
	3000-22	50/CS	XL-RT

Clinical Applications:

- Carpel Tunnel (Overuse Injuries)
- Incomplete Fractures
- Tendonitis
- Soft Tissue Injuries
- Post-Cast Support for Healing



Discover our full range of medical products at medzah.com. Visit us for detailed information and see how our products can enhance your well-being.