

Wrist Brace

Custom Fit | Secure Wrist Support | Functional Hand Movement



Product Overview

The Wrist Brace provides effective wrist support and immobilization while allowing unrestricted movement of the hand and thumb. Its expandable panel conforms to the patient’s anatomy for a customized fit, while the auto-lock strap helps maintain secure positioning. The low-cut palmar design supports comfortable, functional movement during recovery.

Features and Benefits

- Custom Fit: Expandable panel conforms to the patient’s wrist.
- Secure Support: Auto-lock strap maintains the adjusted fit.
- Functional Mobility: Allows unrestricted hand and thumb movement.
- Low-Profile Design: Supports up to 90-degree MP flexion.
- Multiple Options: Available in left- and right-hand configurations.
- Clinical Use: Suitable for carpal tunnel syndrome, incomplete fractures, tendonitis, soft-tissue injuries, and post-cast support.

Product Specifications

- Product Type: Custom-fitting wrist brace
- Closure Type: Adjustable auto-lock strap
- Support Stays: Aluminum and polypropylene
- Main Body Material: Polyester
- Additional Materials: Nylon, steel, and hook-and-loop components
- Laterality: Left and right

Packaging Details

Item No.	QTY	Size
3000-11	50/CS	XS-LT
3000-12	50/CS	XS-RT
3000-10	50/CS	Standard-LT
3000-20	50/CS	Standard-RT
3000-21	50/CS	XL-LT
3000-22	50/CS	XL-RT

