

Thumb Spica

Custom Fit | Secure Thumb Immobilization | Functional Hand Movement



Product Overview

The Thumb Spica is a patented, off-the-shelf support designed to provide a custom fit. Its neutral thumb position allows opposition while fully immobilizing the thumb. The expandable panel automatically conforms to the patient's anatomy, while the auto-lock strap secures the selected fit. The standard size fits virtually all patients, helping reduce the need for excess inventory.

Features and Benefits

- Custom Fit: Expandable panel automatically conforms to the patient's anatomy.
- Secure Immobilization: Neutral thumb positioning allows opposition while fully immobilizing the thumb.
- Auto-Lock Adjustment: Secures the customized sizing and support level.
- Functional Mobility: Low-cut palmar design allows up to 90-degree MP flexion.
- Broad Patient Fit: Standard size fits virtually all patients, helping reduce excess inventory.

Product Specifications

- Product Type: Custom-fitting thumb spica
- Thumb Position: Neutral
- Closure Type: Adjustable auto-lock strap
- Palmar Design: Low-cut at the palmar crease for 90-degree MP flexion
- Main Material: Open-cell polyurethane foam
- Support Stays: Aluminum and nylon
- Additional Materials: Polyester, nylon, steel, and POM
- Laterality: Left and right

Packaging Details

Item No.	QTY	Size
3010-11	50/CS	XS-LT
3010-12	50/CS	XS-RT
3010-10	50/CS	Standard-LT
3010-20	50/CS	Standard-RT
3010-21	50/CS	XL-LT
3010-22	50/CS	XL-RT

