

Air Ankle Stirrup:

Ankle Stirrup Orthosis with Inflatable Air Bladder Liner

Air Ankle Stirrups are designed to provide the stability and relief necessary for a comfortable recovery.

The stirrups feature two adjustable air bladders and low-profile design that can be worn in shoes. Easy-to-use hook and loop closures ensure a quick and hassle-free application.

Key Features:

- Inflatable air bladders provide a comfortable fit and actively reduce edema.
- Hook and loop closures allow for quick and convenient clinical application.
- Bilateral shells provide rigid support, allowing dorsal and plantar flexion while controlling inversion/eversion.

Clinical Applications:

- Acute Ankle Sprains
- Post-Operative Support
- Ankle Injury Rehabilitation
- Chronic Ankle Instability
- Prophylactic Use

MATERIALS USED IN PRODUCT:

| Component | Material | |
|-----------|---|--|
| Sling | Polyester fabric, nylon webbing | |
| Loops | Polypropylene loops | |
| Straps | Nylon webbing, nylon hook and loop w/ plastic rivets | |
| Air Cell | Polyurethane foam Pressure Sensitive adhesive TPU bladder | |
| Tape | Nylon | |
| Stockinet | Polyester | |
| Thread | Polyester | |

AIR ANKLE STIRRUP

| Product | ltem no. | QTY. | Size |
|----------------------|----------|-------|------|
| Air Ankle Stirrup | 4000-00 | 50/CS | 10" |

Discover our full range of medical products at medzah.com. Visit us for detailed information and see how our products can enhance your well-being.