

Air Ankle Stirrups

Ankle Stabilization | Inflatable Air Bladder Liner | Low-Profile Support



Product Overview

Air Ankle Stirrups are designed to provide stability and comfortable support during recovery. Two adjustable air bladders provide a secure, customized fit and help reduce edema. The low-profile design can be worn inside most shoes, while hook-and-loop closures allow quick and convenient clinical application.

Features and Benefits

- Adjustable Air Support: Inflatable air bladders conform comfortably around the ankle.
- Edema Management: Air bladder compression helps reduce swelling during recovery.
- Rigid Stabilization: Bilateral shells provide dependable medial and lateral ankle support.
- Controlled Movement: Allows dorsiflexion and plantar flexion while controlling inversion and eversion.
- Low-Profile Design: Fits comfortably inside most shoes.
- Quick Application: Hook-and-loop closures allow easy fitting, adjustment, and removal.

Product Specifications

- Product Type: Air ankle stirrup orthosis
- Support Design: Bilateral rigid shells
- Liner Type: Two adjustable inflatable air bladders
- Closure Type: Hook-and-loop straps
- Sling Material: Polyester fabric and nylon webbing
- Loop Material: Polypropylene
- Strap Material: Nylon webbing, nylon hook-and-loop material, and plastic rivets
- Air Cell Material: Polyurethane foam, pressure-sensitive adhesive, and TPU bladder
- Tape Material: Nylon
- Stockinette Material: Polyester
- Thread Material: Polyester

Packaging Details

Item No.	QTY	Size
4000-00	50/CS	10 in

