

Round Toe Post Op Shoe

Post-Operative Shoe with Round Toe and Mesh Top

The Round Toe Post-Op Shoe is meticulously designed to offer unparalleled durability and comfort during the recovery phase. This shoe is essential for post-operative care, providing the necessary protection and support to aid in a smooth and effective healing process.

Key Features:

- Extended rocker sole is both rigid and durable.
- Breathable Mesh Upper.
- Reinforced Heel Counter and Support.
- Sole materials reduce sliding.

Clinical Applications:

- Broken Toes
- Metatarsal Fractures
- Foot and Ankle Sprains

MATERIALS USED IN PRODUCT:

Component	Material	
Upper material	Polyester Mesh + PE Foam + Cotton	
Vamp Overlay	PU + Nylon Hook	
Collar Lining	Nylon + PE foam	
Collar Lining Foam	PE Foam	
Heel Counter	TPU	
Vamp Velcro Loop	Polyester	

Round Toe Post Op Shoe

Product	ltem no.	QTY.	Size
Round Toe Post Op Shoe	5000-01	36/CS	S-Male
	5000-02	36/CS	M-Male
	5000-03	36/CS	L-Male
	5000-04	36/CS	XL-Male
	5010-01	36/CS	S-Female
	5010-02	36/CS	M-Female
	5010-03	36/CS	L-Female

Discover our full range of medical products at medzah.com. Visit us for detailed information and see how our products can enhance your well-being.