

Post-Op Shoes

Post-Operative Support | Breathable Comfort | Enhanced Stability



Product Overview

The Post-Op Shoes is designed to provide protection, stability, and comfort throughout the recovery process. Its rigid, durable rocker sole supports stable movement, while the breathable mesh upper and reinforced heel counter help provide a secure and comfortable fit during healing.

Features and Benefits

- Post-Operative Support: Helps protect and support the foot throughout the healing process.
- Rigid Rocker Sole: Provides durable support and promotes stable movement.
- Breathable Upper: Mesh construction promotes airflow and patient comfort.
- Reinforced Heel: Adds stability and provides additional heel support.
- Improved Traction: Sole materials help reduce sliding during walking.
- Adjustable Fit: Hook-and-loop straps allow convenient application and adjustment.

Product Specifications

- Product Type: Round toe post-operative shoe
- Toe Design: Open round toe
- Upper Material: Polyester mesh, PE foam, and cotton
- Vamp Overlay: PU with nylon hook
- Collar Lining: Nylon and PE foam
- Collar Lining Foam: PE foam
- Heel Counter: TPU
- Vamp Loop Material: Polyester
- Closure Type: Hook-and-loop straps
- Available Options: Male and female sizes

Packaging Details

Item No.	QTY	Size
5000-01	36/CS	S-Male
5000-02	36/CS	M-Male
5000-03	36/CS	L-Male
5000-04	36/CS	XL-Male
5010-01	36/CS	S-Female
5010-02	36/CS	M-Female
5010-03	36/CS	L-Female

