

Cervical Collars

Gentle Support | Breathable Comfort | Adjustable Fit



Product Overview

Cervical Collars provide comfortable neck support while helping reduce movement of the head and cervical vertebrae. The open-cell foam construction offers gentle stabilization, while the ribbed stockinet cover promotes air circulation and helps prevent heat buildup and perspiration. An adjustable hook-and-loop closure provides a secure, custom fit during recovery.

Features and Benefits

- Gentle Stabilization: Helps reduce movement of the head and cervical vertebrae.
- Breathable Comfort: Ribbed stockinet cover promotes air circulation and helps prevent heat buildup and perspiration.
- Soft Construction: Open-cell foam provides comfortable support.
- Adjustable Fit: Hook-and-loop closure allows quick and secure adjustment.
- Recovery Support: Suitable for whiplash and neck sprains.

Product Specifications

- Product Type: Cervical collar
- Foam Material: Open-cell polyurethane foam
- Cover Type: Ribbed stockinet
- Cover Material: Polyester-spandex jersey knit fabric
- Closure Type: Hook-and-loop
- Closure Material: Polyester
- Strap Material: Polyester

Packaging Details

Item No.	QTY	Size
6000-25	40/CS	2.5 in
6000-30	40/CS	3 in
6000-40	40/CS	4 in

