

Disposable Surgical Mask

Fluid-Resistant | Medical Grade | Secure Fit



Product Overview

This disposable surgical mask is designed for high-efficiency protection in clinical, hospital, and healthcare settings. Featuring a three-layer construction with a melt-blown filter, it delivers $\geq 98\%$ bacterial filtration and reliable fluid resistance. Comfortable earloop fastenings and hygienic individual packaging make it ideal for daily use in sterile environments.

Features and Benefits

- Medical-Grade Protection: Compliant with Type IIR standards for hospital and clinical use.
- Three-Layer Barrier: Inner and outer non-woven layers with melt-blown filter for advanced filtration.
- Secure Earloop Fit: Ensures comfort and a reliable seal throughout extended wear.
- High Filtration Efficiency: Provides $\geq 98\%$ bacterial filtration with added fluid resistance.
- Hygienic Packaging: Individually wrapped for contamination control and convenience.

Product Specifications

- Type: 3-ply, fluid-resistant (Type IIR)
- Filtration: $\geq 98\%$ BFE
- Compliance: EN 14683:2019, Type IIR
- Packaging: 50/pack, 2,000/case
- Sterilization: Nonsterile

Packaging Details

Item	UOM
BSC-2000	2,000/CS

