

Plastic Specimen Containers

Secure Seal | Tamper-Evident | Lab-Ready



Product Overview

This plastic specimen container features a secure double-click locking seal, providing added confidence in sample collection and transport. Designed to hold up to 90 mL (3 oz), the container is EO (ethylene oxide) sterilized and 95kPa approved, ensuring it meets key clinical safety and transport standards. A tamper-evident seal and temperature strip help safeguard sample integrity, while its durable design ensures reliability across medical, diagnostic, and lab environments.

Features and Benefits

- Double-Click Locking Cap: Provides a secure, leak-proof seal.
- 95kPa Certified: Approved for compliant specimen transport under pressure without risk of leakage.
- EO Sterilized: Ethylene oxide sterilization ensures safe use in clinical and laboratory environments.
- Tamper-Evident Design: Includes integrated seal and temperature strip to confirm sample authenticity.
- Durable Construction: Shatter-resistant design helps prevent cracking and leakage during use.
- Labeled Panel: Pre-printed identification area for easy patient labeling and processing.
- Versatile Use: Ideal for collecting urine, sputum, or other biological specimens across clinical settings.

Product Specifications

- Material: Medical-grade polypropylene
- Capacity: 3 oz / 90 mL
- Seal Type: Double-click locking cap
- Sterilization: EO (Ethylene Oxide)
- Pressure Rating: 95 kPa approved
- Tamper Evident: Yes – with integrated seal and temperature strip
- Label Area: Pre-printed identification panel

Packaging Details

Item No.	Volume	Quantity
EVON101	3OZ / 90ML	400 PCS/CS
EVON102	3OZ / 90ML	400 PCS/CS

