

Aluminium Quad Canes

Height Adjustable | Stable Quad Base | Lightweight Durability



Product Overview

Medzah Aluminum Quad Canes provide sturdy and adjustable support for individuals recovering from injury or managing long-term mobility needs. Constructed from lightweight aluminum and equipped with a stable K-shaped quad base, these canes promote balance, comfort, and confidence while walking.

Features and Benefits

- Quad base design enhances stability, especially on uneven surfaces.
- Shock-absorbing PVC handle provides comfort and reduces hand fatigue.
- Tool-free height adjustment allows for a personalized fit and proper posture.

Product Specifications

- Material: Aluminum
- Adjustable Height: 27–37 inches
- Weight Capacity: 300 lbs
- Tube Wall Thickness: 1.2 mm
- Handle Material: PVC
- Hand Orientation: Left/Right Compatible
- Base Options: Small or Large Quad Base

Packaging Details

Item No.	Size	Weight	Quantity Per Case
MUJI1210-01	Small Base	300 lbs	20
MUJI1210-02	Large Base	300 lbs	10

