

Standard Aluminum Canes

Height Adjustable | Lightweight Design | Foam Comfort Grip



Product Overview

Medzah’s Standard Aluminum Canes are designed for individuals who require reliable support and improved balance in their daily routines. These lightweight canes adjust easily to suit users of various heights and feature a comfortable foam handle with a classic hook-over design for convenience.

Features and Benefits

- Foam grip enhances user comfort and prevents hand fatigue.
- Lightweight aluminum construction ensures easy maneuverability.
- Adjustable to match user height and posture for optimal support.

Product Specifications

- Material: Aluminum
- Adjustable Height: 27–37 inches
- Weight Capacity: 220 lbs
- Tube Wall Thickness: 1.2 mm
- Handle Material: Foam
- Handle Type: Standard (curved)
- Cane Weight: 0.75–0.8 lb

Packaging Details

Item No.	Size	Weight	Quantity Per Case
MUJI3210-0206	Standard	220 lbs	6

