

Perineal Instant Cold Pack

Instant Cold Therapy | Perineal Care | Fold-to-Activate



Product Overview

The Perineal Instant Cold Pack is designed to provide convenient localized cold therapy for perineal care. Its fold-to-activate design breaks the internal water pouch to initiate cooling, allowing the pack to be applied externally to the area of discomfort without prior refrigeration.

Features and Benefits

- Dual Perineal Care: Combines instant cold therapy with absorbent OB-pad functionality for postpartum care.
- Instant Cold Relief: Fold-to-activate design provides convenient, single-use cooling for localized pain and soreness.
- Enhanced Absorbency: Absorbent construction helps manage fluid and reduce leakage during use.
- Simple Application: Easy-to-activate design requires no refrigeration before use.

Product Specifications

- Product Type: Perineal instant cold pack
- Item Number: SAM-4000
- Size: 4.5 x 14.25 in
- Contents: Calcium ammonium nitrate and water
- Activation Method: Fold-to-activate internal water pouch
- Application: External use only
- Use Type: Single-use
- Transport Classification: Non-hazardous for transport

Packaging Details

Item No.	QTY	Size
SAM-4000	24/CS	4.5 x 14.25 in

