

Non-woven Instant Cold

Direct to skin application | Consistent Cooling | Easy Activation



Product Overview

Medzah Non-Woven Instant Cold Packs are designed to provide fast, consistent, and convenient cooling without the need for refrigeration. These single-use packs activate by folding and kneading, delivering approximately 20 minutes of cold therapy while maintaining a safe temperature above 50°F. Made with a durable non-woven pouch, they are ideal for medical, athletic, emergency, and recovery settings where quick cold therapy is needed.

Features and Benefits

- Instant Cooling Relief: Activates on demand to help reduce discomfort from minor injuries, swelling, or soreness.
- No Refrigeration Needed: Stores at room temperature and is ready to use when cold therapy is needed.
- Controlled Temperature: Maintains safe cooling performance without dropping below 50°F.
- Moisture-Free Application: Helps prevent condensation buildup on the skin during use.
- Durable Pouch Design: Non-woven PE exterior supports reliable handling during activation and application.
- Multi-Setting Use: Suitable for clinics, first aid kits, sports programs, schools, and emergency response.

Product Specifications

- Size: 5 x 7.5 inches
- Cooling Duration: Approximately 20 minutes
- Minimum Temperature: Does not drop below 50°F
- Formula Base: Calcium ammonium nitrate and water
- Activation Method: Fold and knead
- Application: Single-use cold therapy

Packaging Details

Item No.	Quantity	Size
SAM-3001	24/CS	5 x 7.5 in

Instructions



Fold product in half



Knead



Apply

