


ADJUSTABLE CRUTCHES:

Adjustable Aluminium Crutches Available in 3 Different Sizes

Medzah’s crutches are expertly designed to provide optimal support for individuals with leg, knee, and ankle injuries.

These crutches aim to alleviate the strain typically associated with walking while injured, allowing users to move with greater ease and comfort. Crafted from high-quality, durable materials, Medzah’s crutches ensure stability and safety, making them an essential aid for recovery and mobility.



MEDZAH'S CRUTCHES: SUPERIOR SUPPORT FOR LEG, KNEE, AND ANKLE RECOVERY

Product	Item no.	Size	Adjustable Height Range	Tube Thickness	Handle Material	Weight	QTY.
Adjustable Aluminium Crutches Available in 3 Different Sizes	LYN-0814-01	Youth	4ft 6in - 5ft 2in	1.1 mm	Rubber	400 lbs	8 Pairs
	LYN-0814-02	Adult	5ft 2in - 5ft 10in	1.2 mm	Rubber	400 lbs	8 Pairs
	LYN-0814-03	Tall	5ft 10in - 6ft 6in	1.2 mm	Rubber	400 lbs	8 Pairs

Explore our full range and experience our full range of medical products on our website. Visit us at www.medzah.com for detailed information and to learn how our products can improve your well-being