

ALCOHOL PREP PADS:Sterile skin prep for procedures.

Where cleanliness meets convenience, ensuring precision in every procedure.

Alcohol prep pads are small, individually-wrapped wipes that are made of non-woven fabric and saturated with 70% isopropyl alcohol. They are used to disinfect the skin prior to medical procedures such as injections, blood draws, or IV insertions. This Alcohol prep pad stands up to standards of USP.



ALCOHOL PREP PADS: The Gold Standard in Sterilization

Product	Item no.	Size	Qty.	Thickness	UOM	Dimensions
ALCOHOL PREP PADS	ALC001	Medium	5000/cs	2 PLY	200/BX 25BX/CS	1x2 in
	ALC002	Large	1000/cs	2 PLY	100/BX 10BX/CS	2x3 in

Explore our full range of incontinence and other medical products on our website. Visit us at www.medzah.com for detailed information and to learn how our products can improve your well-being