

MedWings Brief Diapers

Maximum Absorbency | Adjustable Fit | Odor Control



Product Overview

MedWings Brief Diapers are designed for users who require secure, high-absorbency protection. With a super-absorbent core, soft non-woven topsheet, and breathable backsheet, these briefs deliver comfort and dryness throughout extended wear. Adjustable waist straps provide a snug, custom fit in both standing and lying positions, while built-in odor control adds discretion and confidence in daily use.

Features and Benefits

- Super Absorbent Core: Provides reliable dryness and leak protection.
- Adjustable Waist Straps: Allows a custom fit for both standing and lying application.
- Odor Control Layer: Helps maintain discretion and freshness.
- Soft & Breathable Design: Promotes airflow and reduces skin irritation.
- Secure Fit: Designed for heavy-duty protection in clinical or home settings.

Product Specifications

- Absorbency Level: Maximum
- Top Layer: Soft non-woven comfort layer
- Core: High-capacity absorbent core
- Backsheet: Breathable, skin-safe barrier
- Material: Latex-free
- Use: Disposable, single-use underpad
- Weight Capacity: Up to 300 lbs

Packaging Details

Item No.	Size	Hip Dimensions	Quantity	Weight Capacity
IPAD03	Medium	31.50–41.34 in	96/CS	100-170 lbs
IPAD04	Large	35.43–47.24 in	72/CS	170-240 lbs
IPAD05	X-large	≥43.31 in	60/CS	240-300 lbs

