

Quick-Collect Urine Cup



Reduces Infection



Prevents Medical Accidents



Increases Work Efficiency



Improved Customer Satisfaction



Easy To Transfer in Bulk

Feature & Benefit

- All in one urine test cup with a integrated, easily detachable tube
- Allows collection of multiple samples through tube replacement
- Enhances workplace satisfaction and promotes a safer inspection environment
- Simplifies sample processing with an attached barcode

Competitor Products



Quick-Collect Urine Cup



Quick-Collect Urine Cup

How it works

The **Quick-Collect Urine Cup** simplifies urine collection with an integrated detachable tube. Patients provide a urine sample directly into the cup, which then easily pours into the attached tube. The tube can be detached for immediate transport to the lab, minimizing transfer steps, reducing the risk of spills and contamination, and streamlining the overall testing process.

Features: Safe and Trustworthy, coming with Speed.

Shortened sample processing time

Works worth 4 hours per 1,000 unit can be shortened up to 50 minutes (based on sample processing time)

Safe working environment

Urine is transferred into the airtight container from the urine test cup, preventing contamination of the surroundings and reducing the risk of infection for workers.

Accurate barcode attachment

The barcode containing the patient's information is attached to the tube rather than the cup, preventing medical accidents such as samples getting mixed together.

Better evaluation from the patients

Reducing testing times allows patients to spend less time in the hospital, increasing customer satisfaction (Shortened TAT).

Facilitating Mass Transportation

While existing large-volume cups presented transportation challenges, the smaller tube format enables the efficient transfer of numerous samples

Product Efficiency

The Quick-Collect Urine Cup is an innovative product designed to improve the efficiency and safety of urine sample collection in medical environments.

All-in-One Design

Combines a container and a detachable tube, making it easy to collect and handle urine samples

Safe and Reliable

The urine is transferred into an airtight container, reducing the risk of contamination and infection

Efficient Processing

The attached barcode simplifies sample processing and reduces the risk of sample mix-ups

Mass Transportation

Smaller tubes enhance the efficiency of transporting large quantities of samples

Improved Patient Experience

Patients benefit from shorter hospital stays thanks to faster processing times

Ordering Information

Product	Item no.	Cup Volume	QTY.
Quick-Collect Urine Cup	EVON001	3OZ / 90ML	240 PCS/CS

