

MEDZAH WALKER:

Adjustable Folding Walkers

Our lightweight and durable walkers are designed to provide reliable support for individuals with mobility limitations.

Constructed with either aluminum or steel, these walkers offer adjustable heights to accommodate various user needs. The comfortable grip handles and sturdy frame ensure secure and easy use. The easy lock mechanism allows for quick and convenient opening and closing.



Explore our full range and experience our full range of medical products on our website. Visit us at www.medzah.com for detailed information and to learn how our products can improve your well-being



CONTACT US

(860)-288-8696
 info@medzah.com
 www.medzah.com
 99 Commerce Street
 Glastonbury, CT 06033

ORDERING INFORMATION:

Adjustable Folding Walkers

Adjustable Steel Wheeled Walker

Specification	Description
SKU	MDZH904SE
Adjustable Height Range	32.5 In - 38.5 In
Weight	400 lbs
Quantity Per Case	1
Wheels	5 In
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Steel

Adjustable Steel Wheeled Walker

Specification	Description
SKU	MDZH904S
Adjustable Height Range	32.5 In - 38.5 In
Weight	400 lbs
Quantity Per Case	4
Wheels	5 In
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Steel

Adjustable Steel Walker

Specification	Description
SKU	MDZH804S
Adjustable Height Range	32 In - 38 In
Weight	400 lbs
Quantity Per Case	4
Wheels	N/A
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Steel

Adjustable Aluminum Walker

Specification	Description
SKU	MDZH804
Adjustable Height Range	32 In - 38 In
Weight	300 lbs
Quantity Per Case	4
Wheels	N/A
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Aluminum

Adjustable Aluminum Youth Wheeled Walker

Specification	Description
SKU	MDZH704
Adjustable Height Range	24 In - 30 In
Weight	300 lbs
Wheels	3 In
Quantity Per Case	4
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Aluminum

Adjustable Aluminum Wheeled Walker

Specification	Description
SKU	MDZH904
Adjustable Height Range	32.5 In - 38.5 In
Weight	300 lbs
Quantity Per Case	4
Wheels	5 In
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Aluminum



CONTACT US

(860)-288-8696
 info@medzah.com
 www.medzah.com
 99 Commerce Street
 Glastonbury, CT 06033