



Introducing our precision-engineered elbow support—a cutting-edge solution for swift recovery and enhanced mobility. This innovative immobilization device is expertly crafted to safeguard IV sites while allowing complete freedom of movement. Made with premium nylon-laminated foam and reinforced by robust aluminum supports, it ensures superior durability and comfort. The design is completed with intuitive hook-and-loop closures for effortless application and removal, providing both ease of use and reliable performance. Elevate your recovery experience with a product that combines protection, comfort, and flexibility.

Features and Benefits

- Elbow immobilization solution designed to protect IV sites on the arm
- Allows unrestricted arm movement while ensuring site safety
- Made from nylon-laminated foam for comfort
- Reinforced with durable aluminum supports for stability
- Easy-to-use hook-and-loop tab closures for quick application and removal

Product	ltem no.	Size	Length	Width	Unit
ELBOW SPLINT	Med-1001	Extra Small	7″	4"	Each
	Med-1002	Small	9″	6"	Each
	Med-1003	Medium	11.5″	7″	Each
	Med-1004	Large	15.5″	10"	Each

