

Elbow Splint

IV-Site Protection | Comfortable Support | Easy Application



Product Overview

The Elbow Splint is designed to immobilize the elbow and protect IV sites on the arm while allowing unrestricted arm movement. Constructed from comfortable nylon-laminated foam and reinforced with durable aluminum supports, it provides dependable stability and support. Hook-and-loop tab closures allow quick application, adjustment, and removal.

Features and Benefits

- IV-Site Protection: Helps safeguard IV sites located on the arm.
- Unrestricted Movement: Allows arm movement while helping keep the IV site protected.
- Comfortable Construction: Nylon-laminated foam provides soft and comfortable support.
- Reliable Stability: Durable aluminum supports reinforce the splint and help maintain its shape.
- Quick Application: Hook-and-loop tab closures allow easy fitting, adjustment, and removal.
- Multiple Size Options: Available in four sizes to accommodate different patient requirements.

Product Specifications

- Product Type: Elbow immobilization splint
- Intended Use: Protection of IV sites on the arm
- Primary Material: Nylon-laminated foam
- Reinforcement Material: Aluminum supports
- Closure Type: Hook-and-loop tab closures
- Movement: Allows unrestricted arm movement

Packaging Details

Item No.	Size	Length	Width	Unit
MED-1001	Extra Small	7 in	4 in	Each
MED-1002	Small	9 in	6 in	Each
MED-1003	Medium	11.5 in	7 in	Each
MED-1004	Large	15.5 in	10 in	Each

