

Offset Handle Fashion Canes

Height Adjustable | Lightweight Support | Stylish Design



Product Overview

Medzah Offset Handle Fashion Canes are thoughtfully crafted for users who value both function and style. Made from lightweight yet durable aluminum, these canes adjust from 27 to 37 inches to suit a range of users and support up to 220 lbs. With colorful patterns and a soft grip, they offer a comfortable and modern alternative to traditional chrome walking aids.

Features and Benefits

- Stylish Appeal: Unique designs provide a fashionable alternative to institutional-looking aids.
- Comfortable Grip: Handle ensures all-day comfort and includes wrist strap for convenience.
- Adjustable Fit: Height range accommodates a variety of user preferences and postures.

Product Specifications

- Material: Aluminum
- Adjustable Height: 27–37 inches
- Weight Capacity: 220 lbs
- Tube Wall Thickness: 1.2 mm
- Handle Material: Polypropylene (PP)
- Cane Weight: 0.8 lb (0.4 kg)
- Size: Standard

Packaging Details

Item No.	Size	Weight	Quantity Per Case
MUJI3210-0106	Standard	220 lbs	6

