

Instant Hot Pack

Consistent Heat | Portable Relief | Easy Activation



Product Overview

Medzah Instant Hot Packs are designed to provide fast, consistent, and controlled therapeutic heat without the need for external heating. These single-use packs activate with a simple squeeze and deliver approximately 20 minutes of soothing warmth between 105°F and 110°F. Built for convenient heat therapy, they are ideal for home care, clinical, athletic, and recovery settings where quick, reliable warmth is needed. Use with a protective layer or towel between the pack and skin for safe application.

Features and Benefits

- Consistent Heat Output: Delivers controlled therapeutic warmth between 105°F and 110°F.
- Easy Activation: Simple fold and knead design activates heat quickly when needed.
- 20-Minute Duration: Provides reliable heat therapy for short treatment sessions.
- Safe Formulation: Sodium acetate and water-based design supports consistent performance
- Long Shelf Life: 3-year shelf life helps support inventory planning and storage.

Product Specifications


- Size: 6 x 6 inches
- Heat Duration: 20 minutes
- Temperature Range: 105°F to 110°F
- Shelf Life: 3 years

Packaging Details


Item No.	Quantity	Size
SAM-2000	24/CS	6 x 6 in

Medzah's Hot Pack Vs Current Markets Hot Pack

MEDZAH'S Instant Hot Pack Sodium Acetate Test



Average temperature during single heat therapy treatment³



Instructions



Fold product in half



Knead



Apply

