

Non-woven Instant Hot Pack

Consistent Heat | Direct-to-Skin Comfort | Easy Activation



Product Overview

Medzah Non-Woven Instant Hot Packs are designed to provide fast, consistent, and controlled therapeutic heat without the need for external heating. These single-use packs activate with a simple squeeze and deliver approximately 20 minutes of soothing warmth between 105°F and 110°F. Made with a soft, skin-friendly non-woven exterior, they are comfortable for direct-to-skin application and ideal for home care, clinical, athletic, or recovery settings where convenient heat therapy is needed.

Features and Benefits

- Consistent Heat Output: Delivers controlled therapeutic warmth between 105°F and 110°F.
- Direct-to-Skin Application: Soft non-woven pouch allows comfortable use without an extra cover.
- Easy Activation: Simple fold and knead design activates heat quickly when needed.
- 20-Minute Duration: Provides reliable heat therapy for short treatment sessions.
- Safe Formulation: Sodium acetate and water-based design supports consistent performance
- Long Shelf Life: 3-year shelf life helps support inventory planning and storage.

Product Specifications

- Size: 6 x 6 inches
- Heat Duration: 20 minutes
- Temperature Range: 105°F to 110°F
- Shelf Life: 3 years

Packaging Details

Item No.	Quantity	Size
SAM-2001	24/CS	6 x 6 in

Instructions



Fold product in half




Knead



Apply

Medzah's Hot Pack Vs Current Markets Hot Pack

MEDZAH'S Instant Hot Pack Sodium Acetate Test



Average temperature during single heat therapy treatment³

