

Instant Cold Packs

Instant Cooling | Portable Relief | First Aid Ready



Product Overview

Medzah Instant Cold Packs are designed to provide fast, convenient cold therapy without the need for refrigeration. These single-use packs activate by folding and kneading, creating an endothermic cooling reaction that delivers approximately 20 minutes of relief while maintaining a controlled temperature above 50°F. Compact and easy to store, they are ideal for clinics, first aid kits, athletic programs, schools, workplaces, and emergency response settings.

Features and Benefits

- Calcium ammonium nitrate & water base for consistent controlled cooling
- CNo moisture accumulation on the skin
- Duration for 20 minutes
- Never drops below 50°F

Product Specifications

- Size: 4.75 x 7 inches
- Cooling Duration: Approximately 20 minutes
- Minimum Temperature: Does not drop below 50°F
- Formula Base: Calcium ammonium nitrate and water
- Activation Method: Fold and knead
- Application: Single-use cold therapy

Packaging Details

Item No.	Quantity	Size
SAM-3000	16/CS	4.75 x 7 in

Instructions



Fold product in half



Knead



Apply

