

# INFANT HEEL WARMER:

## Instant heat therapy

**Our Infant Heel Warmer provides a safe and effective means of providing localized warmth to infant feet.**

The device utilizes a food-grade sodium acetate and water mixture that generates a consistent temperature between 101°F (38°C) and 104°F (40°C) upon activation. This temperature range is ideal for promoting comfort and aiding in blood flow, making it particularly useful for procedures requiring heel sticks or for soothing discomfort associated with cold environments.

**How it works:** The unit is activated by squeezing firmly on the inner fluid pouch. Rapid crystallization occurs when the liquid contents are exposed to the minute crystals of sodium acetate contained within the outer pouch. The exothermic reaction causes the unit to heat up.

Reaches up to  
104° F



5.5"

3.75"

### INFANT HEEL WARMER

Product	Item no.	QTY.	Size
INFANT HEEL WARMER	SAM-108	25ea/bx 4bx/cs	5.5"x3.75"

Explore with us and experience our full range of medical products on our website. Visit us at [www.medzah.com](http://www.medzah.com) for detailed information and to learn how our products can improve your well-being