

Splash Mask Face Shields

UltraClear Vision | Anti-Fog Design | Latex-Free Safety



Product Overview

The Splash Mask Face Shield is engineered to deliver unmatched optical clarity and full-face fluid protection in clinical and procedural environments. Featuring dual-sided anti-fog, anti-static, and anti-glare coatings, these latex-free shields ensure clear visibility and safe use even during high-intensity procedures. Lightweight and comfortable, each shield fits securely with either elastic or Velcro straps and is available in vented or solid headpiece options.

Features and Benefits

- Anti-Fog Coating: Prevents visual obstruction in humid or extended use environments.
- Anti-Static Design: Reduces static buildup for added safety.
- Latex-Free Materials: Safe for users with latex sensitivities.
- Crystal-Clear Shield: Maintains sharp visibility for precise work.
- Anti-Glare Surface: Minimizes visual distortion under bright lights.

Product Specifications

- Shield Type: Full or half coverage
- Headpiece Options: Vented or solid
- Strap Type: Elastic or Velcro
- Coatings: Anti-fog, anti-static, anti-glare (both sides)
- Material: Latex-free, fiberglass-free
- Filtration: ≥99% @ 0.1 micron
- Packaging: 2 bags/case × 25 units/bag (50 per case)

Packaging Details

Item No.	Shield	Headpiece	Strap	Quantity/Cs
MEDIT-000	Full	Vented	Velcro	2 bags/case x 25 units/bag
MEDIT-001	Full	Vented	Elastic	2 bags/case x 25 units/bag
MEDIT-002	Full	Solid	Velcro	2 bags/case x 25 units/bag
MEDIT-003	Full	Solid	Elastic	2 bags/case x 25 units/bag
MEDIT-005	Half	Vented	Velcro	2 bags/case x 25 units/bag
MEDIT-006	Half	Vented	Elastic	2 bags/case x 25 units/bag



Contact us today at: Email: info@medzah.com | Phone: (860) 288-8696 | Website: medzah.com

MDMEDIT-000 | 071025