

# Stretch Vinyl Gloves

Eco-Friendly | Chemotherapy-Tested | Clinical-Grade Quality



## Product Overview

These disposable vinyl stretch gloves offer a perfect blend of comfort, protection, and a ordability for medical/food handling tasks. Made from lightweight vinyl, they provide a close,flexible fit that allows for natural dexterity and tactile sensitivity. Ideal for medical, food and general purpose tasks. These single-use gloves promote hygiene and simplify cleanup.



## Features and Benefits

- Thin and Flexible: Comfortable fit for natural hand movement.
- Medical Grade Glove: Tested in accordance with ASTM F1671 & ASTM D5250
- Disposable: Minimizes contamination risk and simplifies cleanup.
- Food Safe: Meets relevant standards for safe food handling.
- Stretchy Comfort: Universal fit for easy on/off.
- Cost-Effective: Ideal for various food handling tasks

## Product Specifications: Physical Properties:

Finger Thickness: 4.3mil	Tensile Strength (MPA): min. 11
Palm Thickness: 3.5mil	Elongation (%): min. 300
Length: 9.5in	Property: ASTM D5250

## Sizing & Packaging Details

Reference No.	Description	Quantity
VLMD01	Powder Free Vinyl Stretch-XS	100/BX 10BX/CS
VLMD02	Powder Free Vinyl Stretch-S	100/BX 10BX/CS
VLMD03	Powder Free Vinyl Stretch-M	100/BX 10BX/CS
VLMD04	Powder Free Vinyl Stretch-L	100/BX 10BX/CS
VLMD05	Powder Free Vinyl Stretch-XL	100/BX 10BX/CS