

PRODUCT CATALOG

2026





10 Million +
products sold

Nationwide
shipping support
across the U.S.

99%+
Customer
Satisfaction Rate

A Customer Centric Solutions Oriented Medical Supply Company

Our commitment is to deliver the highest quality medical supplies to our customers. We recognize that each customer has unique needs, and we create solutions to offer specialized services to meet those needs.



Excellent Quality



More Than 10
Locations Worldwide



Over 150 Highly Motivated,
Knowledgeable Employees



MBE Certified



Gnemsdc Awarded Fastest
Growing Business of 2024



Awarded Greater New England's
Manufacturer of the year for 2025

U.S.-Based Support

A responsive North American team that makes communication, quoting, sampling, and issue resolution easier for healthcare customers.

Custom Kitting & Packaging

Experience building tailored medical kits, private label packaging, and product-specific solutions for labs, clinics, and care programs.

Global Manufacturing Network

Access to international manufacturing and sourcing capabilities, helping customers balance quality, cost, and scalability.

Flexible, Customer-Focused Service

Smaller and more hands-on than large distributors, allowing Medzah to move faster, support custom requests, and adapt to customer needs.

A Product for every **CARE** need



Examination Gloves

Biodel Nitrile Gloves
Bionex Nitrile Gloves
MediServe Nitrile Gloves
Nite Shield Nitrile Gloves
Stretch Vinyl Gloves



Mobility Aids

Walkers
Crutches
Canes



Facial Protection

Splash Mask Face Shields
Clear Shieldz
Protective Glasses
Protective Goggles
Disposable 3-Ply Face Masks



Incontinence

MedEase Pullup Diapers
Medmax Underpads
MedWings Brief Diapers



Diagnostics

Blood Culture Collection Kits
Oncology Kits
Urine Collection Containers
Plastic Specimen Containers



Pain Relief

Refillable Ice Bags with Ties
Non-Woven Instant Cold Packs
Instant Cold Packs
Non-Woven Instant Hot Packs
Instant Hot Packs
Infant Heel Warmers



Orthopedic Soft Goods

Elbow Splints
Air Ankle Stirrups
Wrist Braces
Post Op Shoes
Compression Knee Immobilizers
Shoulder Immobilizers
Thumb Spicas
Cervical Collars
Tri-Panel Knee Immobilizer



Wound Care

Gauze Pads
Alcohol Prep Pads



Clinical Packaging

Cold Gel Pack
Insulated Box Liners
Corrugated Shipping Boxes

Comprehensive Medical Supply Solutions

for Every Care Setting



Examination Gloves



Biodel Nitrile Gloves
Bionex Nitrile Gloves
MediServe Nitrile Gloves
Nite Shield Nitrile Gloves
Stretch Vinyl Gloves



Mobility Aids



Walkers
Crutches
Canes



Facial Protection



Splash Mask Face Shields
Clear Shieldz
Protective Glasses
Protective Goggles
Disposable 3-Ply Face Masks



Incontinence



MedEase Pullup Diapers
Medmax Underpads
MedWings Brief Diapers



Specimen Container



Urine Collection Containers
Plastic Specimen Containers



Pain Relief



Refillable Ice Bags with Ties
Non-Woven Instant Cold Packs
Instant Cold Packs
Non-Woven Instant Hot Packs
Instant Hot Packs
Infant Heel Warmers



Orthopedic Soft Goods



Elbow Splints
Air Ankle Stirrups
Wrist Braces
Post Op Shoes
Compression Knee Immobilizers
Shoulder Immobilizers
Thumb Spicas
Cervical Collars
Tri-Panel Knee Immobilizer



Wound Care



Biodel Nitrile Gloves
Bionex Nitrile Gloves
MediServe Nitrile Gloves
Nite Shield Nitrile Gloves
Stretch Vinyl Gloves

Personal Protection

Medzah Examination Gloves are designed to support reliable hand protection across clinical, laboratory, and general healthcare environments. This category includes nitrile and vinyl glove options suited for routine exams, patient care, infection control, and everyday medical use, offering a range of comfort, durability, and protection levels to meet different facility needs.



PERSONAL PROTECTION

Biodel Nitrile Gloves

Specification	Description
Color	Violet blue
Material	Nitrile
Sizes	XS-XL
Packaging	200/bx 10bx/cs
Grade	Medical Grade
Finger Thickness	4.0 mil
Palm Thickness	3.0 mil
Texture	Fingertip Textured
Latex	Latex Free
Powder	Powder Free
Cuff Length	9.5 in

Item No.	Description	Pkg
BDMD01	Biodel Gloves (XS)	200/bx 10bx/cs
BDMD02	Biodel Gloves (S)	200/bx 10bx/cs
BDMD03	Biodel Gloves (M)	200/bx 10bx/cs
BDMD04	Biodel Gloves (L)	200/bx 10bx/cs
BDMD05	Biodel Gloves (XL)	200/bx 10bx/cs

Nite Shield Nitrile Gloves

Specification	Description
Color	Black
Material	Nitrile
Sizes	S-XXL
Packaging	100/bx 10bx/cs
Grade	Medical Grade
Finger Thickness	5.5 mil
Palm Thickness	3.9 mil
Texture	Fingertip Textured
Latex	Latex Free
Powder	Powder Free
Cuff Length	9.5 in

Item No.	Description	Pkg
NSMD02	Niteshield Gloves (S)	100/bx 10bx/cs
NSMD03	Niteshield Gloves (M)	100/bx 10bx/cs
NSMD04	Niteshield Gloves (L)	100/bx 10bx/cs
NSMD05	Niteshield Gloves (XL)	100/bx 10bx/cs
NSMD06	Niteshield Gloves (XXL)	90/bx 10bx/cs

Stretch Vinyl Gloves

Specification	Description
Color	Yellow
Material	Vinyl
Sizes	XS-XL
Packaging	100/bx 10bx/cs
Grade	Medical Grade
Finger Thickness	4.3 mil
Palm Thickness	3.5 mil
Texture	Fingertip Smooth
Latex	Latex Free
Powder	Powder Free
Cuff Length	9.5 in

Item No.	Description	Pkg
VLMD01	Vinyl Gloves (XS)	100/bx 10bx/cs
VLMD02	Vinyl Gloves (S)	100/bx 10bx/cs
VLMD03	Vinyl Gloves (M)	100/bx 10bx/cs
VLMD04	Vinyl Gloves (L)	100/bx 10bx/cs
VLMD05	Vinyl Gloves (XL)	100/bx 10bx/cs



PERSONAL PROTECTION

Bionex Nitrile Gloves

Specification	Description
Color	Violet blue
Material	Nitrile
Sizes	XS-XL
Packaging	200/bx 10bx/cs (XL- 180/bx 10bx/cs)
Grade	Medical Grade
Finger Thickness	4.0 mil
Palm Thickness	3.0 mil
Texture	Fingertip Textured
Latex	Latex Free
Powder	Powder Free
Cuff Length	9.5 in

Item No.	Description	Pkg
BNMD01	Bionex Gloves (XS)	200/bx 10bx/cs
BNMD02	Bionex Gloves (S)	200/bx 10bx/cs
BNMD03	Bionex Gloves (M)	200/bx 10bx/cs
BNMD04	Bionex Gloves (L)	200/bx 10bx/cs
BNMD05	Bionex Gloves (XL)	180/bx 10bx/cs

MediServe Nitrile Gloves

Specification	Description
Color	Violet Blue
Material	Nitrile
Sizes	XS-XL
Packaging	250/bx 10bx/cs
Grade	Medical Grade
Finger Thickness	3.5 mil
Palm Thickness	2.5 mil
Texture	Fingertip Textured
Latex	Latex Free
Powder	Powder Free
Cuff Length	9.5 in

Item No.	Description	Pkg
MSMD01	MediServe Gloves (XS)	250/bx 10bx/cs
MSMD02	MediServe Gloves (S)	250/bx 10bx/cs
MSMD03	MediServe Gloves (M)	250/bx 10bx/cs
MSMD04	MediServe Gloves (L)	250/bx 10bx/cs
MSMD05	MediServe Gloves (XL)	230/bx 10bx/cs



Mobility Aids

Medzah Mobility Aids are designed to support safe movement, stability, and patient independence. This category includes walkers, crutches, and canes for rehabilitation, recovery, long-term mobility support, and everyday patient assistance.



MOBILITY AIDS

Product #MDZH904SE

Specifiction	Description
Adjustable Height Range	32.5 In - 38.5 In
Weight	400 lbs
Quantity Per Case	1
Wheels	5 In
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Steel

Product #MDZH904S


Specifiction	Description
Adjustable Height Range	32.5 In - 38.5 In
Weight	400 lbs
Quantity Per Case	4
Wheels	5 In
Tube Thickness	1.2 mm
Handle Material	PVC
Material	Steel

Product #MDZH704

Specifiction	Description
Adjustable Height Range	24 In - 30 In
Weight	300 lbs
Wheels	3 In
Quantity Per Case	4
Tube Thickness	1.2 mm
Handle Material	Rubber
Material	Aluminum

Adjustable Aluminium Crutches Available in 3 Different Sizes

Item no.	Size	Adjustable Height Range	Tube Thickness	Handle Material	Weight	QTY.
LYN-0814-01	Youth	4ft 6in - 5ft 2in	1.1 mm	Rubber	400 lbs	8 Pairs
LYN-0814-02	Adult	5ft 2in - 5ft 10in	1.2 mm	Rubber	400 lbs	8 Pairs
LYN-0814-03	Tall	5ft 10in - 6ft 6in	1.2 mm	Rubber	400 lbs	8 Pairs



MOBILITY AIDS

ADJUSTABLE STANDARD CANE

Product #MUJI3210-0206

Specification	Description
Size	Standard
Adjustable Height Range	27-37 in
Weight	220 lbs
Quantity Per Case	6
Tube Thickness	1.2 mm
Handle Material	Foam



ADJUSTABLE OFFSET HANDLE CANE

Product #MUJI3210-0106

Specification	Description
Size	Standard
Adjustable Height Range	27-37 in
Weight	220 lbs
Quantity Per Case	6
Tube Thickness	1.2 mm
Handle Material	PP



ADJUSTABLE QUAD CANE

Product #MUJI1210-01

Specification	Description
Size	Small Base
Adjustable Height Range	27-37 in
Weight	300 lbs
Quantity Per Case	20
Tube Thickness	1.2 mm
Handle Material	PVC



Product #MUJI1210-02

Specification	Description
Size	Large Base
Adjustable Height Range	27-37 in
Weight	300 lbs
Quantity Per Case	10
Tube Thickness	1.2 mm
Handle Material	PVC



Facial Protection

Medzah Facial Protection products are designed to help reduce exposure to splashes, droplets, and airborne particles in medical and general care settings. This category includes face shields, protective eyewear, goggles, and disposable masks for practical frontline protection.



FACIAL PROTECTION

Clear Sheildz

Item no.	Description	QTY.	Sterilization
MEDIT-500	Pre Assembled	100/CS	Nonsterile
MEDIT-506	Dispenser Box Assemble Clear Sheildz	120/CS	Nonsterile
MEDIT-501	Eye Protection Lens 100 EA/BX, 5BX/CS	500/CS	Nonsterile
MEDIT-502	Eye Protection Frames 25EA/BX, 8BX/CS	200/CS	Nonsterile

- Wrap around frame for comfort, a wrap around design for protection
- Top down protection PATENTED self closing holes seal when attached
- Finger print free handling with the colored top strip
- Sealed lens mounting for added fluid protection
- Clip-on/off design for easy changing
- Natural rubber-free for universal use
- Anti fog, anti glare and anti static



Splash Mask Face Shields

Item no.	Shield	Headpiece	Strap	Sterilization
MEDIT-000	Full	Vented	Velcro	2 bags/case x 25 units/bag
MEDIT-001	Full	Vented	Elastic	2 bags/case x 25 units/bag
MEDIT-002	Full	Solid	Velcro	2 bags/case x 25 units/bag
MEDIT-003	Full	Solid	Elastic	2 bags/case x 25 units/bag
MEDIT-005	Half	Vented	Velcro	2 bags/case x 25 units/bag
MEDIT-006	Half	Vented	Elastic	2 bags/case x 25 units/bag

- **Anti-fog coating (both sides):** Ensures clear vision in any environment
- **Anti-static:** Reduces the buildup of static electricity for added safety
- **Latex-safe:** Thoughtfully designed for those with latex sensitivities
- **Optically clear:** Provides outstanding clarity during use
- **Anti-glare:** Minimizes glare for enhanced visibility



FACIAL PROTECTION

Protective Glasses

Item no.	Description	QTY.	Sterilization
MEDIT-300	Protective Glasses	100/CS	Nonsterile

- Accommodates prescription glasses: Designed for a comfortable fit over your eyewear.
- 99.9% UV protection: Provides superior safety for your eyes against harmful UV rays.
- Side ventilation: Offers added comfort and minimize fogging for clear vision.
- High-quality optical lenses: Ensures exceptional clarity during use.
- Latex-safe: Crafted with care for those with latex sensitivities.
- Natural rubber-free for universal use
- Scratch resistant.



Protective Goggles

Item no.	Description	QTY.	Sterilization
MEDIT-328	Protective Goggles	100/CS	Nonsterile


- Adjustable elastic strap: Easily customize the fit for optimal comfort and protection.
- 99.9% UV protection: Offers superior defense for your eyes against harmful UV rays
- High-quality optical lenses: Ensures outstanding clarity during use.
- Latex Safe: Mindfully crafter for those with latex sensitivities.



Disposable 3-Ply Face Masks

Item no.	Product	UOM
MDZH-FM01	Disposable 3-Ply Face Masks	2,000/CS

- Medical Use: Designed for clinical, hospital, and healthcare settings requiring ASTM Level 1 protection.
- Construction: Premium 3-layer construction using 25gsm non-woven and melt-blown layers.
- Filtration: Provides $\geq 95\%$ BFE and $\geq 95\%$ PFE.
- Fastenings: Mindfully crafter for those with latex sensitivities.
- Material Safety: Latex-free and hypoallergenic for sensitive skin.
- Fluid Resistance: Synthetic blood resistance up to 80 mmHg.
- Breathability: < 5.0 mm H₂O/cm² for comfortable wear.
- Flammability: Class 1.
- Sterilization: Non-sterile.
- Type: 3-ply, ASTM Level 1 (Low fluids).
- Packaging: 2,000/CS (50/BX, 40BX/CS).



Incontinence

Medzah Incontinence products are designed to support patient comfort, hygiene, and everyday care needs. This category includes pull-up diapers, brief diapers, and underpads for use in healthcare facilities, long-term care, and home care settings.




INCONTINENCE PRODUCTS

MedWings Brief Diapers

Item no.	Size	Hip Dimensions	QTY	Weight Capacity
IPAD03	Medium	31.50–41.34 in.	96/cs	100-170 lbs
IPAD04	Large	35.43–47.24 in.	72/cs	170-240 lbs
IPAD05	X-large	≥43.31 in.	60/cs	240-300 lbs

- High absorbency adult wing diaper
- Superior sap for lasting dryness
- Plush material for softness
- Advanced odor control
- Breathable layers



MedEase PullUp Diapers

Item no.	Size	Waist Dimension	QTY	Weight Capacity
IPPD03	Medium	24 x 45 in.	80/cs	100-170 lbs
IPPD04	Large	27 x 53 in.	72/cs	170-240 lbs
IPPD05	X-large	32 x 61 in.	56/cs	240-300 lbs

- Soft back sheet for comfort
- Secure fit and discreetly designed
- Easy off sides
- High absorbency core



Medmax Underpads

Item no.	Size	QTY	Weight Capacity
IPUP02	23 x 36 in	150/cs	N/A
IPUP03	30 x 36 in.	40/cs	400lbs
IPUP04	26 x 36 in.	40/cs	N/A
IPUP05	30 x 36 in.	40/cs	N/A

- Cutting-edge design with super-absorbent core for lasting dryness
- Soft non-woven topsheet for enhanced comfort
- Breathable backsheet to promote airflow and skin health



Specimen Container

Medzah Specimen Container support organized specimen collection, testing workflows, and healthcare procedure preparation. This category includes collection kits, oncology kits, urine containers, and specimen containers designed for reliable clinical and laboratory use.



SPECIMEN CONTAINER

Plastic Specimen Container

Specification	Description
Item no.	EVON101
Volume	3OZ / 90ML
QTY.	400 PCS/CS

- EO sterilization ensures a high level of cleanliness and suitability for medical and lab use
- Tamper-evident seal to guarantee sample integrity and prevent unauthorized access
- Generous 90 ml capacity to accommodate an ample amount of specimens
- Single click lock mechanism for a secure and leak-proof seal, 95kPa



Urine Collection Container

Specification	Description
Item no.	CP10221
Cup Volume	8 OZ Cup
Biodegradable	No
QTY.	500 Sets

- Ensures that samples remain in tact during transport
- Made of high quality paper at 280 gsm
- Container is sturdy and leak proof
- 18 gsm plastic liner for cup
- 2 gsm plastic lid



Pain Relief

Medzah Pain Relief products are designed to provide convenient hot and cold therapy for recovery, injury care, and patient comfort. This category includes instant hot packs, cold packs, ice bags, and infant heel warmers for clinical and home care use.



PAIN RELIEF

Instant Cold Pack

Item no.	Description	QTY	Size
SAM-3000	Instant Cold Pack	16/cs	4.75x7 In

- **Contents:** Calcium ammonium nitrate and water, kept separate until activation
- **Activation Method:** Squeeze the pack to break the internal barrier
- **Chemical Reaction:** Calcium ammonium nitrate dissolves in water, triggering an endothermic reaction
- **Heat Absorption:** The reaction absorbs heat from the surroundings, causing a rapid temperature drop
- **Optimal Cooling Temperature:** Reached quickly
- **Duration:** Provides about 20 minutes of cooling



Instant Hot Pack

Item no.	Description	QTY	Size
SAM-2000	Instant Hot Pack	24/cs	6x6 In

- **Consistent, Therapeutic Warmth:** Provides warmth for up to 20 minutes
- **Unique Blend:** Utilizes sodium acetate and water
- **Safe Temperature Range:** Maintains 105°F to 110°F
- **Activation Method:** Fold the pouch to activate the inner water package
- **Rapid Chemical Reaction:** Generates soothing heat quickly
- **Reliable and Effective:** Offers relief from discomfort without exceeding 110°F



PAIN RELIEF

Non Woven Instant Cold Pack

Item no.	Description	QTY	Size
SAM-3001	Non Woven Instant Cold Pack	24/cs	5x7.5 In

- **Non-Woven Fabric:** Ensures even distribution of the cooling effect and offers comfortable direct application to the skin
- **Contents:** Contains calcium ammonium nitrate and water, separated until activation
- **Activation Method:** Squeeze the pack to break the internal barrier, allowing the calcium ammonium nitrate to dissolve in the water
- **Endothermic Reaction:** Absorbs heat from the surroundings, resulting in a rapid drop in temperature
- **Optimal Cooling Temperature:** Reached almost instantly
- **Duration:** Provides approximately 20 minutes of cooling
- **Ideal Uses:** Medical use, sports injuries, and first aid kits

Non-woven Instant Hot Pack

Item no.	Description	QTY	Size
SAM-2001	Non-woven Instant Hot Pack	24/cs	6x6 In

- **Non-Woven Instant Hot Pack:** Designed with your comfort in mind
- **Skin-Friendly Exterior:** Allows for direct application
- **Quick and Effective Heat Therapy**
- **Safe Temperature Range:** 105°F to 110°F
- **Duration:** Provides soothing warmth for up to 20 minutes
- **Relieves Discomfort Effortlessly**



PAIN RELIEF

Instant Cold Pack

Item no.	Description	QTY	Size
SAM-1001	Instant Cold Pack	24/CS	6×9 In

- Innovative metal disk mechanism guarantees instant fluid release
- Convenient fold over and squish activation for hassle-free use
- Ammonium nitrate technology for fast and reliable cooling
- Compact and portable design, perfect for on-the-go relief



PAIN RELIEF

Infant Heel Warmer

Item no.	Description	QTY	Size
SAM-108	Infant Heel Warmer	25ea/bx 4bx/cs	5.5x3.75 In

- **Generates consistent heat:** Between 101°F (38°C) and 104°F (40°C)
- **Activation Method:** Squeeze the inner fluid pouch to initiate an exothermic reaction
- **Aids in comfort**
- **Promotes blood flow**
- **Useful for heel stick procedures**
- **Soothes discomfort in cold environments**



Refillable Ice Bags With Ties

Item no.	Description	QTY	Size
SAM-5001	Refillable Ice Bag, Large with Ties	50/CS	5.5"x3.75"

- **Large 6.5" x 14" Size with Secure Ties:** Stays securely in place during cold therapy.
- **Refillable & Reusable:** Clear design for easy ice-level monitoring and quick refills.
- **Single-Patient Use** Helps maintain hygiene in clinical and home care settings.
- **Latex-Free** Safe for users with latex sensitivities.
- **25-Count Pack** Ideal for hospitals, clinics, and home care.



Orthopedic Soft Goods

Medzah Orthopedic Soft Goods provide supportive bracing and immobilization solutions for injury recovery, post-operative care, and joint stabilization. This category includes braces, splints, collars, immobilizers, and orthopedic supports for a range of patient needs.



Orthopedic Soft Goods

Elbow Splints

Item no.	Size	Length	Width	Unit
MED-1001	Extra Small	7"	4"	Each
MED-1002	Small	9"	6"	Each
MED-1003	Medium	11.5"	7"	Each
MED-1004	Large	15.5"	10"	Each

- Elbow immobilization solution designed to protect IV sites on the arm
- Allows unrestricted arm movement while ensuring site safety
- Made from nylon-laminated foam for comfort
- Reinforced with durable aluminum supports for stability
- Easy-to-use hook-and-loop tab closures for quick application and removal



Air Ankle Stirrups

Item no.	QTY	Size
4000-00	50/CS	10"

- Inflatable air bladders provide a comfortable fit and actively reduce edema.
- Hook and loop closures allow for quick and convenient clinical application.
- Bilateral shells provide rigid support, allowing dorsal and plantar flexion while controlling inversion/eversion.



Wrist Braces

Item no.	QTY	Size
3000-11	50/CS	XS-LT
3000-12	50/CS	XS-RT
3000-10	50/CS	Standard-LT
3000-20	50/CS	Standard-RT
3000-21	50/CS	XL-LT
3000-22	50/CS	XL-RT

- Standard size fits virtually all patients, reducing need for excess inventory.
- Panel automatically expands, providing a custom off-the-shelf fit.
- Auto-lock strap secures custom sizing.
- Cut low at the palmar crease to allow for 90-degree MP flexion.



Orthopedic Soft Goods

Post Op Shoes

Item no.	QTY	Size
5000-01	36/CS	S-Male
5000-02	36/CS	M-Male
5000-03	36/CS	L-Male
5000-04	36/CS	XL-Male
5010-01	36/CS	S-Female
5010-02	36/CS	M-Female
5010-03	36/CS	L-Female

- Extended rocker sole is both rigid and durable.
- Breathable Mesh Upper.
- Reinforced Heel Counter and Support.
- Sole materials reduce sliding.



Compression Knee Immobilizers

Item no.	QTY	Size
2000-00	12/CS	12"
2000-01	12/CS	16"
2000-02	12/CS	20"
2000-03	12/CS	24"

- Contoured, reinforced stays with center "I" beam provide rigidity and a custom fit.
- Soft, foam laminate material adds comfort and resists staining. Can be easily trimmed to size.
- Adjustable elastic straps for even compression.
- Universal design fits most.



Shoulder Immobilizers

Item no.	QTY	Size
1000-00	40/CS	Extra-Small
1000-01	40/CS	Small
1000-02	40/CS	Medium
1000-03	40/CS	Large
1000-04	40/CS	Extra-Large

- Durable cotton sling ensures long-lasting use and optimal support.
- Hook and loop closure adjustments on all straps allow for a personalized fit.
- Fits both the right and left arm, adding to its versatility.
- The waist strap helps prevent shoulder rotation, promoting effective healing.
- Included shoulder pad offers additional comfort, enhancing patient compliance.



Orthopedic Soft Goods

Thumb Spicas

Item no.	QTY	Size
3010-11	50/CS	XS-LT
3010-12	50/CS	XS-RT
3010-10	50/CS	Standard-LT
3010-20	50/CS	Standard-RT
3010-21	50/CS	XL-LT
3010-22	50/CS	XL-RT

- Standard size fits virtually all patients, reducing need for excess inventory.
- Panel automatically expands, providing a custom off-the-shelf fit.
- Auto-lock strap secures custom sizing.
- Cut low at the palmar crease to allow for 90-degree MP flexion.

Cervical Collars

Item no.	QTY	Size
6000-25	40/CS	2.5"
6000-30	40/CS	3"
6000-40	40/CS	4"

- Offers comfortable support while reducing head and cervical vertebrae movement.
- Ribbed stockinet cover promotes air circulation and helps prevent heat buildup and perspiration.
- Hook and loop adjustment provides a custom fit.

Tri-Panel Knee Immobilizer

Item no.	QTY	Size
40012	40/CS	Extra-Small
40016	40/CS	Small
40020	40/CS	Medium
40024	40/CS	Large

- Effective support for knee immobilization during recovery
- Comfortable fit to prevent chafing and irritation
- Adjustable design to accommodate leg size and swelling changes
- Easy to put on and take off, even for users with limited mobility
- Durable construction for long-term use without losing support
- Comfortable fit to prevent chafing and irritation
- Adjustable design to accommodate leg size and swelling changes
- Easy to put on and take off, even for users with limited mobility
- Durable construction for long-term use without losing support



Wound Care

Medzah Wound Care products provide dependable support for basic wound management, preparation, and patient care. This category includes gauze pads and alcohol prep pads designed for everyday clinical, first aid, and healthcare facility use.



WOUND CARE

Item no.	Size	Qty.	Thickness	UOM	Dimensions
ALC001	Medium	5000/cs	2 PLY	200/BX 25BX/CS	1x2 in
ALC002	Large	1000/cs	2 PLY	100/BX 10BX/CS	2x3 in

- **Sterile:** individually wrapped to ensure safety and free of contaminants
- **Effective:** 70% isopropyl alcohol kills bacteria and microorganisms
- **Affordable:** cost-effective for healthcare facilities
- **Fast-acting:** disinfects skin in just a few seconds
- **Convenient:** small, portable, and easy to use



Gauze Pads			
Item no.	Description	Thickness	QTY.
GCL111	Woven non sterile gauze 2x2	12 PLY	8000 PCS
GCL114	Woven non sterile gauze 2x2	8 PLY	5000 PCS
GCL115	Woven sterile gauze 2x2	8 PLY	3000 PCS

- Versatile for different medical applications and sizes
- Non-linting to prevent debris or fibers in the wound
- Highly absorbent for effective wound care
- Soft and comfortable for patient use



NEW

CLINICAL PACKAGING SOLUTIONS

Medzah Clinical Packaging solutions are designed to help maintain product integrity during storage and shipment. This category includes cold gel packs, insulated box liners, corrugated shipping boxes, and custom packaging solutions for temperature-sensitive medical and laboratory products.



Learn more at medzah.com



CONTACT US

860-288-8696

info@medzah.com

medzah.com

99 Commerce Street
Glastonbury, CT 06033

SCAN TO LEARN MORE

