


# Clear Shieldz


## Use in three easy steps



1

### Position the Shieldz


Place the Clear Shieldz over the eyes and position the wrap-around frame comfortably.



2

### Secure the Fit

Adjust the Clear Shieldz for a secure fit and ensure the self-closing top seal is in place.



3

### Check Visibility

Confirm the lens is clear and provides unobstructed visibility before use.

### Reusable Design:

Repeated uses. Applies to Medzah Clear Shieldz, Item Numbers: MEDIT-500, MEDIT-506.

### Contents:

Clear Shieldz eyewear with frame and lens, assembled. Supplied in a dispenser box. Nonsterile.

### Caution:

Inspect the Clear Shieldz before use. Do not use if damaged or if visibility is impaired. Ensure a secure fit before use.

### Storage and Handling:

Store the Clear Shieldz in a clean, dry location to maintain lens clarity and product condition. Handle the shieldz carefully to prevent scratches, smudges, or other damage to the lens. Avoid placing heavy objects on the Clear Shieldz.

### Important Use Information:

Position the Clear Shieldz securely over the eyes and ensure the wrap-around frame provides a comfortable fit. The shieldz is compatible with prescription glasses and provides fluid protection through its patented self-closing top seal. Handle the lens carefully when removing or changing it.

