


Gauze Pads

Application in three easy steps



1

Prepare the Gauze


Remove the gauze pad from its packaging and prepare it for use.



2

Apply to Wound

Place the gauze pad gently over the affected area for wound protection and absorption.



3

Secure the Dressing

Secure the gauze with tape or a bandage as directed by a healthcare professional.

Single-Use Design:

Single use. Applies to Medzah Gauze Pads, Item Number: GCL111 (Non-Sterile, 12-ply), GCL114 (8-ply), GCL115 (Sterile, 8-ply).

Contents:

Woven gauze made from soft, pliable synthetic woven material.

Caution:

Do not use if the package or gauze is damaged or contaminated. Use with appropriate wound-care procedures and follow a healthcare professional's directions. Do not use if the sterile package is opened or damaged.

Storage and Handling:

Store in a clean, dry location at room temperature. Keep packaged gauze protected from moisture, contamination, and damage. Inspect the packaging before use.

Important Use Information:

Inspect the gauze and packaging before use. Do not use if damaged or contaminated. Apply gently over the wound to absorb exudate and provide protection. Follow a healthcare professional's directions.

