


MedEase Pull-Up Diapers


Application in two easy steps



1

Position the Diaper

Step into the diaper and pull it up around the waist with the absorbent core aligned correctly.



2

Adjust the Fit

Adjust the diaper around the waist and hips to ensure a secure, comfortable fit.

Single-Use Design:

Single use. Applies to Medzah MedEase Pull-Up Diapers, Item No. IPPD03 (M), IPPD04 (L), IPPD05 (XL).

Contents:

Soft, high-absorbency core with a stretchable elastic waistband and breathable design.

Caution:

Do not reuse the diaper. Ensure it is properly fitted around the waist and hips to help prevent leakage, maintain comfort, and avoid excessive pressure on the skin.

Storage and Handling:

Store the diapers in a clean, dry location at room temperature. Keep them protected from excessive moisture, heat, direct sunlight, and sharp objects that may damage the packaging or product.

Important Use Information:

Select the size by waist and hip measurements. Pull the diaper up and adjust it around the waist and hips for a secure fit, and ensure the diaper provides reliable protection without causing excessive pressure or discomfort.

