


MedWings Brief Diapers


Application in three easy steps



1

Position the Diaper


Position the diaper around the waist with the absorbent core aligned correctly.



2

Adjust the Fit

Adjust the diaper around the waist and hips to ensure a secure, comfortable fit.



3

Secure Side Panels

Secure and adjust side panels to provide a snug fit and help prevent leakage.

Single-Use Design:

Single use. Applies to Medzah MedWings Brief Diapers, Item Numbers: IPAD03 (M), IPAD04 (L), IPAD05 (XL).

Contents:

Super-absorbent core with soft non-woven topsheet, breathable backsheet, and adjustable waist straps.

Caution:

For external use only. Do not reuse. Do not use if the package or brief is damaged. Discontinue use if skin irritation occurs.

Storage and Handling:

Store MedWings Brief Diapers in a clean, dry location at room temperature. Keep away from excessive moisture, heat, direct sunlight, and sharp objects that may damage the product or packaging.

Important Use Information:

Select the appropriate size according to hip dimensions and weight capacity. Position the brief correctly, adjust the side panels for a secure fit, and ensure the brief remains comfortably positioned during use.

