


Nite Shield Nitrile Gloves


Application in three easy steps



1

Select the Right Size


Choose a glove size that fits your hand snugly and comfortably.



2

Don the Glove

Gently pull the glove on and fit it smoothly over the fingers, palm, and wrist.



3

Remove and Dispose

Carefully peel the glove off inside out and discard it after a single use.

Single-Use Design:

Single use. Applies to Medzah Nite Shield Nitrile Gloves, Item Numbers: NSMD02 (S), NSMD03 (M), NSMD04 (L), NSMD05 (XL), NSMD06 (XXL).

Contents:

Black nitrile examination glove, powder-free and latex-free. Length 9.5 in; finger thickness 0.13 mm (5.5 mil); palm thickness 0.10 mm (4.0 mil).

Caution:

Single use only. Do not reuse. Do not use if the glove or package is damaged. Discontinue use if skin irritation occurs.

Storage and Handling:

Store Nite Shield Nitrile Gloves in a clean, dry location at room temperature. Keep away from excessive moisture, heat, direct sunlight, and sharp objects that may damage the product or packaging.

Important Use Information:

Select the appropriate size for a snug, comfortable fit. Tested for use with chemotherapy drugs per ASTM D6978-05; refer to the spec sheet for the drugs tested. Dispose of used gloves according to local regulations.

