



Protective Glasses

Use in three easy steps



1 Position the Glasses

Place the protective glasses over your eyes and prescription eyewear for full-wrap coverage.



2 Adjust the Fit

Ensure the glasses sit comfortably and securely over your existing eyewear.



3 Check Visibility

Confirm the lenses provide clear, undistorted vision and the side vents remain unobstructed.

Reusable Design:

Repeated uses. Applies to Medzah Protective Glasses, Item Number: MEDIT-300.

Contents:

Pre-assembled Clear Shieldz (frame with lens). Nonsterile.

Caution:

Inspect the glasses before each use to ensure the lenses and frame are in good condition. Do not use if the product is damaged or visibility is impaired. Ensure the glasses fit securely before use.

Storage and Handling:

Store in a clean, dry location when not in use. Handle the glasses carefully to avoid scratching or damaging the lenses. Keep the ventilation openings clear.

Important Use Information:

These glasses are designed to fit over most prescription eyewear while providing full-wrap protection. Ensure the optical-grade lenses remain clear and the side ventilation is unobstructed during use.

