


Protective Goggles


Use in three easy steps



1

Position the Goggles

Place the goggles securely over your eyes and ensure they sit comfortably.



2

Adjust the Strap

Adjust the elastic strap for a secure, customized, comfortable fit.



3

Check the Fit

Ensure the goggles are properly positioned and the lenses provide clear visibility.

Reusable Design:

Repeated uses. Applies to Medzah Protective Goggles, Item Number: MEDIT-328.

Contents:

Protective goggles with clear optical-grade lens and adjustable elastic strap. 99.9% UV protection. Latex-free. Nonsterile.

Caution:

Inspect the goggles before each use to ensure the lenses, frame, and elastic strap are in good condition. Do not use if any part is damaged or if visibility is impaired. Ensure the goggles are securely fitted before use.

Storage and Handling:

Store the goggles in a clean, dry location when not in use. Handle them carefully to prevent scratches or damage to the optical lenses. Keep the elastic strap and lenses protected from unnecessary wear.

Important Use Information:

Adjust the elastic strap to achieve a secure and comfortable fit around the eyes. Ensure the goggles remain properly positioned during use while maintaining clear visibility through the optical-grade lenses.

